CC DISK Zr
CC DISK Zr HT

CC DISK Zr HT Multilayer
CC DISK Zr SMILE Multilayer

(SL) Navodila za uporabo / (EN) Directions for use / (DE) Gebrauchsanweisung / (CS) Návod k použití
(HR) Uputa za uporabu / (ES) Instrucciones de uso / (RU) Инструкции по использованию


Proizvajalec / Producer / Hersteller
Interdent d.o.o. · SI - 3000 CELJE
Opekarniška cesta 26 · T: +386 (0)3 425-62-00
E: info@interdent.cc · www.interdent.cc

CE 0197

Verzija: 11/2023
Datum: 11.07.2023

Made in Slovenia

Indikacije

	Luska	Inlay & onlay	Reducirana krona	Polno anatomska krona	Reducirana prevleka	Polno anatomski anteriorni most (3 členi)	Polno anatomski posteriorni most (3 členi)	Polno anatomski most (≤7 členov)	Polno anatomski most (≤14 členov)	Abutment
* CC DISK Zr Bela opačna										
CC DISK Zr HT Predbarvana & Bela										
CC DISK Zr HT Multilayer										
CC DISK Zr SMILE Multilayer										

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Indication guide

	Veneer	Inlay & onlay	Reduced crown	Full contour crown	Copyng
* CC DISK Zr White opaque					
CC DISK Zr HT Preshade & White					
CC DISK Zr HT Multilayer					
CC DISK Zr SMILE Multilayer					

* Because of its lower translucency, the CC DISK Zr is suitable for frameworks which will be layered with ceramic material.

Full contour anterior bridge 3(unit)	Full contour posterior bridge 3(unit)	Full contour bridge (≤7unit)	Full contour bridge (≤14unit)	Abutment

CC DISK Zr (white, opaque), CC DISK Zr HT (pure and preshade, translucent), CC DISK Zr HT Multilayer (multilayer preshade, high translucent), CC DISK Zr SMILE Multilayer (multilayer preshade, super translucent) are discs made out of biocompatible pre-sintered zirconium dioxide stabilized with yttrium oxide. They are used in CAD/CAM milling machines. CC DISK Zr (HT, HT Multilayer) comply with demands of the standard EN ISO 6872, type II, class 5. CC DISK Zr Smile Multilayer comply with demands of the standard EN ISO 6872, type II, class 4.

Contraindications: Using in case of bruxism, inadequate strength of tooth structure or insufficient space is not recommended. Usage is also not recommended in case of inadequate oral hygiene or inadequate coordination of masticatory movements. Prosthetic substitutes made out of discs must not be used in case of unsatisfactory preparation of teeth and in case of mistakes in prosthetic substitute.

Intended users: Product is intended to be used by professionals – dental technicians.

Patient selection criteria: Totally or partially edentulous patients.

Physical and chemical characteristics	CC DISK Zr	CC DISK Zr HT	CC DISK Zr HT preshade	CC DISK Zr HT Multilayer	CC DISK Zr SMILE Multilayer
Colour	White opaque	White translucent	high translucent	Multilayer preshade: high translucent	Multilayer preshade: super translucent
Sintered density [g/cm ³]	>6,02	>6,02	>6,02	>6,02	> 6,02
Flexural strength [MPa]	1400±100	1250±100	1250±100	900-1100	600-900
Thermal expansion [10 ⁻⁶ K ⁻¹]	10,5	10,5	10,5	10,5	10,5
Translucency [%]	42	43	43	46	49
Radioactivity [Bq/g]	< 0,1	< 0,1	< 0,1	< 0,1	< 0,1
Solubility [µg/cm ²]	< 50	< 50	< 50	< 50	< 50

Demands for preparation: Doctor of dental medicine is responsible for proper tooth preparation. Design of frame for ceramic veneering must follow the anatomical shape of the tooth. The connectors must be round in shape. Please consider that the connector section has to be dimensioned larger, depending on the construction.

Below written minimal wall thickness dimension are reference values:

Posterior section	Occlusal 0,9 mm	Circular 0,5 mm
Anterior section	Incisal/occlusal 0,7mm	Circular 0,5 mm
Abutment crowns of bridge frameworks with 2 pontics	Incisal/occlusal 1,0mm	Circular 0,7 mm

The connector surface should be at least:

Anterior bridge	With one pontic	7 mm ²
Anterior bridge	With two pontics	9 mm ²
Posterior bridge	With one pontic	9 mm ²
Posterior bridge	With two pontics	12 mm ²

Milling: Milling can be done in any dental CAD/CAM machine. Milling machine's manufacturer must inform user about suitable cutters and revolutions per minute for milling CC DISK Zr. Discs must not come into contact with liquid or oil during milling in order to avoid strength loose or staining.

Wet milling of presintered discs is possible but construction must be dried out completely before sintering process. Drying can be done in drying chamber on 80 °C or under infrared light. As an alternative low-temperature firing program can also be selected, which does not harm the construction. Notches, sharp edges or slight damages, also in inner part, can reduce sintered construction strength, therefore it is important that rounded and smooth surface is ensured after milling.

Colouring: If required, the frame may be coloured to desired colour with standard coloring solutions that are suitable for zirconium oxide ceramics. In order to achieve uniform result, it must be ensured that the framework is completely clean, greaseless and dry before coloring. For the correct working procedure, instructions of the coloring solutions manufacturer should be followed.

Sintering

Recommended sintering program

Standard program:

CC DISK Zr (white, HT, HT preshade, HT Multilayer, Zr SMILE Multilayer)

Fast program:

CC DISK Zr (beli, HT, HT preshade)

Units	1-7	above 7	1-3
Heating rate	3°C per minute up to 300°C	3°C per minute up to 300°C	6°C per minute up to 300°C
Maintenance time	without	without	without
Heating rate	5°C per minute up to 1520°C	4 °C per minute up to 1520°C	10 °C per minute up to 1520°C
Final temperature	1520°C	1520°C	1520°C
Maintenance time at the final temperature	90 minutes	90 minutes	60 minutes
Cooling	controlled 4°C / min to 800°C, then cooling in a closed furnace	controlled 3°C / min to 800°C, then cooling in a closed furnace	controlled 18°C / min to 800°C, then cooling in a closed furnace

Fast program is intended to be used for single crowns and bridges with maximum 3 pontics. In order to achieve the maximum predicted transparency of the material, it is recommended to use only a standard program.

Sintering can be done in all standard furnaces that are suitable for zircon oxide ceramic. Follow the instruction for use of the furnace.

Enlargement factor depends on the lot and it is written on the edge of the disc.

We recommend sintering construction in sintering crucible, placed in suitable granular material. One of the suitable granular type are ZrO₂ sintering spheres. If you use crucible with granular material, place the construction with the occlusal side down on the granular material in the crucible. Prevent contacts with other frames or sintering carriers. If you use sintering carriers or wires than make sure that they do not come into contact with the walls of the caps.

Dental constructions and crowns made out of zirconium oxide can be sintered only in suitable furnaces. If other ceramic materials are fired or alloy is sintered in the same furnace, this can influence the sintering behavior or can cause local discoloration.

Ceramic veneering: You can use any ceramic suitable for veneering zirconium oxide ceramic. Follow the instruction for use of ceramic for veneering.

Finishing after final sintering: Any finishing after final sintering must be done with water spray to prevent local overheating or damage the material.

Storing: If discs are stored in original packaging at room temperature than changes in characteristics of the material are not expected, also if stored for a longer period. Do not expose discs to direct sunlight and keep dry.

Disposal of the residual material: Residual part of discs and sintered part of discs are chemically inert. Dispose them in accordance with local regulation.

Safety precaution! Powder is harmful to health. Use suction unit. At manual treatment use protective mask FFP1 and wear goggles.

Notice: Any serious incident, that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

Side effects: Consider possibility of allergic hypersensitivity to the powder.

Warranty! Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can only be considered as standard values. Our product is subject to further development. Therefore, alteration in construction and composition of working instructions is reserved.

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